

## Product datasheet for **HP227448**

### SPNS2 Human qPCR Primer Pair (NM\_001124758)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	124976
Forward Sequence:	TTACTGGCTCCAGCGTGAAGCA
Reverse Sequence:	GGGACCAGAATGAGGATGAGTG
ACCN:	<u>NM_001124758, NM_001124758.1, NM_001124758.2, BC041772, BC065221, NM_001124758.3</u>
UniProt ID:	<u>Q8IVW8</u>
Synonyms:	DFNB115; SLC62A2; SLC63A2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)