

Product datasheet for **HP227456**

ACOT7 Human qPCR Primer Pair (NM_181864)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11332
Forward Sequence:	CTACACCTCCAAGCACTCTGTG
Reverse Sequence:	CCTTGTCACATTCTTCAGCGAC
ACCN:	<u>NM_181864</u> , <u>NM_181864.1</u> , <u>NM_181864.2</u> , <u>BC017365</u> , <u>BC017365.2</u> , <u>BG332214</u> , <u>BM807217</u> , <u>BQ067921</u> , <u>BT006888</u> , <u>BU594651</u> , <u>NM_181864.3</u>
UniProt ID:	<u>O00154</u>
Synonyms:	ACH1; ACT; BACH; CTE-II; hBACH; LACH; LACH1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)