

Product datasheet for **HP227484**

MTMR11 Human qPCR Primer Pair (NM_181873)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10903
Forward Sequence:	GAACACTGTCCAGATTCCTGGC
Reverse Sequence:	CCAAGGACAGAGCTGATTCAGG
ACCN:	NM_181873 , NM_181873.1 , NM_181873.2 , NM_181873.3 , BC032655 , BC064372 , BC094756 , BC100818 , BC100819 , BC100820 , BC100821 , BC166616 , BQ446006
UniProt ID:	A4FU01
Synonyms:	CRA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)