

Product datasheet for **HP227537**

CYFIP1 Human qPCR Primer Pair (NM_014608)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23191
Forward Sequence:	TACGAGACGCTGCTGAAGCAGA
Reverse Sequence:	TCGTCCAATCGCCAGTTCTAGG
ACCN:	<u>NM_014608</u> , <u>NM_014608.1</u> , <u>NM_014608.2</u> , <u>NM_014608.4</u> , <u>BC005097</u> , <u>BC001306</u> , <u>BC063869</u> , <u>BF979687</u>
UniProt ID:	<u>Q7L576</u>
Synonyms:	P140SRA-1; SHYC; SRA-1; SRA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)