

Product datasheet for **HP227563**

Histone acetyltransferase MYST3 (KAT6A) Human qPCR Primer Pair (NM_006766)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7994
Forward Sequence:	CTTCAGTGAGAGCAGCGAGGAG
Reverse Sequence:	GTGGTGTTTGCCTTTTCGGACT
ACCN:	<u>NM_006766</u> , <u>NM_006766.1</u> , <u>NM_006766.2</u> , <u>NM_006766.3</u> , <u>NM_006766.4</u> , <u>BC018011</u> , <u>BC142659</u> , <u>BC142959</u> , <u>BC172379</u> , <u>BP248857</u> , <u>BP271264</u> , <u>BP348860</u> , <u>NM_006766.5</u>
UniProt ID:	<u>Q92794</u>
Synonyms:	ARTHS; MOZ; MRD32; MYST-3; MYST3; RUNXBP2; ZC2HC6A; ZNF220
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)