

Product datasheet for **HP227643**

PPM1B Human qPCR Primer Pair (NM_002706)

Product data:

| | |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 5495 |
| Forward Sequence: | GGCTTGAGGTATCTGATGACCTG |
| Reverse Sequence: | CACCGCTTCATCTGAGACCTTG |
| ACCN: | NM_002706 , NM_002706.1 , NM_002706.2 , NM_002706.3 , NM_002706.4 , NM_002706.5 , BC064381 , BC064381.1 , BC012002 , BF979498 , BQ773881 , NM_002706.6 |
| UniProt ID: | O75688 |
| Synonyms: | PP2C-beta; PP2C-beta-X; PP2CB; PP2CBETA; PPC2BETAX |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)