

Product datasheet for **HP227648**

ANO6 Human qPCR Primer Pair (NM_001142680)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	196527
Forward Sequence:	TACCTGACTCCACAGACAGCCA
Reverse Sequence:	ATCAGTCTGGGTCCTGGGAGT
ACCN:	NM_001142680 , NM_001142680.1 , BC098410 , BC098410.1 , BC042383 , BC063576 , BC073891 , BC136445
UniProt ID:	Q4KMQ2
Synonyms:	MGC104751; TMEM16F
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)