

Product datasheet for **HP227653**

DDX19B Human qPCR Primer Pair (NM_007242)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11269
Forward Sequence:	GCTTCTCCAAGGAGTCTATGCC
Reverse Sequence:	CAGTACCAGACTGAGATTGGGC
ACCN:	NM_007242 , NM_007242.1 , NM_007242.2 , NM_007242.3 , NM_007242.4 , NM_007242.5 , BC003626 , BC010008 , BC010008.2 , BC110805 , BI756088 , NM_007242.7
UniProt ID:	Q9UMR2
Synonyms:	DBP5; DDX19; RNAh
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)