

Product datasheet for **HP227662**

Profilin 2 (PFN2) Human qPCR Primer Pair (NM_053024)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5217
Forward Sequence:	TCGGCTACTGCGACGCCAAATA
Reverse Sequence:	GATCACTGAGCATTCTTCGCGC
ACCN:	NM_053024 , NM_053024.1 , NM_053024.2 , NM_053024.3 , BC043646 , BC043646.1 , BC002964 , BC018049 , BC095444 , BX463350 , NM_053024.4
UniProt ID:	P35080
Synonyms:	D3S1319E; PFL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)