

Product datasheet for **HP227750**

Paxillin (PXN) Human qPCR Primer Pair (NM_002859)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5829
Forward Sequence:	CTGATGGCTTCGCTGTCCGATT
Reverse Sequence:	GCTTGTTCAAGTCAGACTGCAG
ACCN:	NM_002859 , NM_002859.1 , NM_002859.2 , NM_002859.3 , BC003647 , BC006840 , BC010498 , BC015217 , BC052611 , BC136787 , BC136794 , BC144410 , BM742662 , BU500419 , BX648777
UniProt ID:	P49023
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)