

Product datasheet for **HP227796**

CNTF Receptor alpha (CNTFR) Human qPCR Primer Pair (NM_147164)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1271
Forward Sequence:	GCTTCTACTGCAGCTGGCATCT
Reverse Sequence:	CGAATGTGGCAGCGTTCTTGA
ACCN:	<u>NM_147164</u> , <u>NM_147164.1</u> , <u>NM_147164.2</u> , <u>BC001492</u> , <u>BC001492.2</u> , <u>BI545014</u> , <u>BI829871</u> , <u>BM714356</u> , <u>BT019824</u> , <u>BX364434</u> , <u>NM_147164.3</u>
UniProt ID:	<u>P26992</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)