

Product datasheet for **HP227803**

RYK Human qPCR Primer Pair (NM_002958)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6259
Forward Sequence:	CAGCAAGACCTGGTACACATGG
Reverse Sequence:	CAAGTCTCTGGAGAGGGCATTG
ACCN:	<u>NM_002958</u> , <u>NM_002958.1</u> , <u>NM_002958.2</u> , <u>NM_002958.3</u> , <u>BC021700</u> , <u>BC148614</u> , <u>BC153090</u> , <u>BQ003586</u> , <u>NM_002958.4</u>
UniProt ID:	<u>P34925</u>
Synonyms:	D3S3195; JTK5; JTK5A; RYK1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)