

Product datasheet for **HP227807**

PAK4 Human qPCR Primer Pair (NM_005884)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10298
Forward Sequence:	GATGATTCGGGACAACCTGCCA
Reverse Sequence:	AGGAATGGGTGCTTCAGCAGCT
ACCN:	NM_005884 , NM_005884.1 , NM_005884.2 , NM_005884.3 , BC011368 , BC011368.1 , BC025282 , BC025282.1 , BC002921 , BC034511 , BG773480 , BM129363 , NM_005884.4
UniProt ID:	O96013
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)