

Product datasheet for **HP227825**

FANCL Human qPCR Primer Pair (NM_018062)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55120
Forward Sequence:	GGAGTGCAACAGCACGCAGAAT
Reverse Sequence:	CTGCTCAGCTTAATTCACAGGG
ACCN:	NM_018062 , NM_018062.1 , NM_018062.2 , NM_018062.3 , BC054517 , BC054517.1 , BC009042 , BC037570 , BC041627 , BQ574618 , NM_018062.4
UniProt ID:	Q9NW38
Synonyms:	FAAP43; PHF9; POG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)