

Product datasheet for **HP227837**

WDR4 Human qPCR Primer Pair (NM_033661)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10785
Forward Sequence:	CGTCTACTCCTTTTCGGTGCTG
Reverse Sequence:	GCAGTGAGGATGAAGCGGTCAT
ACCN:	NM_033661 , NM_033661.1 , NM_033661.2 , NM_033661.3 , NM_033661.4 , BC006341 , BC006341.2 , BC001074 , BC018866 , BG722541 , BQ447230 , BU541429 , BU558089
UniProt ID:	P57081
Synonyms:	GAMOS6; hWH; MIGSB; TRM82; TRMT82; Wuho
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)