

Product datasheet for **HP227892**

AP3D1 Human qPCR Primer Pair (NM_003938)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8943
Forward Sequence:	TCGTCACAGAGGAGATGCCTGA
Reverse Sequence:	GCAGTTTCTCGCTGTCGGCTAA
ACCN:	NM_003938 , NM_003938.1 , NM_003938.2 , NM_003938.3 , NM_003938.4 , NM_003938.5 , NM_003938.6 , BC005142 , BC010065 , BC022837 , BM982261 , BQ227203 , NM_003938.8
UniProt ID:	O14617
Synonyms:	ADTD; hBLVR; HPS10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)