

Product datasheet for **HP227899**

Ataxin 3 (ATXN3) Human qPCR Primer Pair (NM_004993)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4287
Forward Sequence:	TCGGAAGAGACGAGAAGCCTAC
Reverse Sequence:	AAGTGCTCCTGAACTGGTGGCT
ACCN:	NM_004993 , NM_004993.1 , NM_004993.2 , NM_004993.3 , NM_004993.4 , NM_004993.5 , BC022245 , BC033711 , BC095402 , BU190081 , BU663113 , NM_004993.6
UniProt ID:	P54252
Synonyms:	AT3; ATX3; JOS; MJD; MJD1; SCA3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)