

Product datasheet for **HP228015**

RNF185 Human qPCR Primer Pair (NM_001135824)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	91445
Forward Sequence:	ACACAGCCAAGGATGCCGTCAT
Reverse Sequence:	GGATGACCTTGTCTCGGCTGAT
ACCN:	<u>NM_001135824</u> , <u>NM_001135824.1</u> , <u>BC035684</u> , <u>BC035684.1</u> , <u>BC012817</u> , <u>BC009504</u> , <u>BC017302</u> , <u>BC026040</u> , <u>BC033166</u> , <u>BX648019</u>
UniProt ID:	<u>Q96GF1</u>
Synonyms:	BSK65-MONO1; BSK65-MONO2; BSK65-PANC1; BSK65-PANC2; BSK65-TEST1; BSK65-TEST2; BSK65-TEST3; FLJ38628; ring finger protein 185
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)