

Product datasheet for **HP228022**

TLE3 Human qPCR Primer Pair (NM_005078)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7090
Forward Sequence:	CCACCATGAACTCGATCACAGAG
Reverse Sequence:	CTTGGCTTCCATGCTGTAGTCC
ACCN:	BC041831 , NM_005078 , NM_005078.1 , NM_005078.2 , NM_005078.3 , BC041831.1 , BC012004 , BC015421 , BC015729 , BC043247 , BM676381 , BP373048 , NM_005078.4
UniProt ID:	Q04726
Synonyms:	ESG; ESG3; GRG3; HsT18976
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)