

## Product datasheet for **HP228023**

### KCNMB2 Human qPCR Primer Pair (NM\_181361)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10242
Forward Sequence:	GGACTGGCTATGATGGTGTGCT
Reverse Sequence:	TTCCGTGATGGACGATTACGC
ACCN:	<u><a href="#">NM_181361</a></u> , <u><a href="#">NM_181361.1</a></u> , <u><a href="#">NM_181361.2</a></u> , <u><a href="#">BC017825</a></u> , <u><a href="#">BC017825.1</a></u> , <u><a href="#">BF434939</a></u> , <u><a href="#">BG701449</a></u> , <u><a href="#">BM983559</a></u> , <u><a href="#">NM_181361.3</a></u>
UniProt ID:	<u><a href="#">Q9Y691</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)