

Product datasheet for **HP228107**

hnRNP L (HNRNPL) Human qPCR Primer Pair (NM_001533)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3191
Forward Sequence:	GTGTGGTGAAGCAGACCTTGT
Reverse Sequence:	CAAACCCACAGTGCTTGCTC
ACCN:	<u>NM_001533</u> , <u>NM_001533.1</u> , <u>NM_001533.2</u> , <u>BC069184</u> , <u>BI256815</u> , <u>BI461733</u> , <u>BM997815</u> , <u>BQ334900</u>
UniProt ID:	<u>P14866</u>
Synonyms:	hnRNP-L; HNRPL; P/OKcl.14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)