

Product datasheet for **HP228220**

LAF4 (AFF3) Human qPCR Primer Pair (NM_002285)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3899
Forward Sequence:	CACCTGGCAAAGTGAACCAAC
Reverse Sequence:	CGATGTGCCATTGTCTGGACTC
ACCN:	<u>NM_002285</u> , <u>NM_002285.1</u> , <u>NM_002285.2</u> , <u>BC021283</u> , <u>BC036895</u> , <u>BC136579</u> , <u>BC144266</u> , <u>BI906088</u> , <u>BU740375</u> , <u>NM_002285.3</u>
UniProt ID:	<u>P51826</u>
Synonyms:	KINS; LAF4; MLLT2-like
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)