

Product datasheet for **HP228223**

CDC27 Human qPCR Primer Pair (NM_001256)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	996
Forward Sequence:	ACACCTCCTGTAATTGATGTGCC
Reverse Sequence:	GGAGTTACCTCTCGGCTATTTCC
ACCN:	<u>NM_001256</u> , <u>NM_001256.1</u> , <u>NM_001256.2</u> , <u>NM_001256.3</u> , <u>BC011656</u> , <u>BI548077</u> , <u>BP234964</u> , <u>NM_001256.6</u>
UniProt ID:	<u>P30260</u>
Synonyms:	ANAPC3; APC3; CDC27Hs; D0S1430E; D17S978E; H-NUC; HNUC; NUC2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)