

## Product datasheet for **HP228241**

### HCCS Human qPCR Primer Pair (NM\_005333)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 3052   |
| Forward Sequence: | GTACGTGGAGTGTCCCATTAGG   |
| Reverse Sequence: | CTGCTCTCGGAATGGATGACTC   |
| ACCN:             | <a href="#">NM_005333</a> , <a href="#">NM_005333.1</a> , <a href="#">NM_005333.2</a> , <a href="#">NM_005333.3</a> , <a href="#">NM_005333.4</a> , <a href="#">BC001691</a> , <a href="#">BC095455</a> , <a href="#">NM_005333.5</a>  |
| UniProt ID:       | <a href="#">P53701</a>   |
| Synonyms:         | CCHL; LSDMCA1; MCOPS7; MLS   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)