

Product datasheet for **HP228281**

VPS53 Human qPCR Primer Pair (NM_001128159)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55275
Forward Sequence:	GCATGGTGCAATGCTCTCAGCT
Reverse Sequence:	TCAGTCCTCCACTGCTGGTTGT
ACCN:	<u>NM_001128159</u> , <u>NM_001128159.1</u> , <u>NM_001128159.2</u> , <u>BC006116</u> , <u>BC006116.1</u> , <u>BC029560</u> , <u>BC039432</u> , <u>BC040223</u> , <u>NM_001128159.3</u>
UniProt ID:	<u>Q5VIR6</u>
Synonyms:	HCCS1; hVps53L; PCH2E; pp13624
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)