

Product datasheet for **HP228307**

LYRM1 Human qPCR Primer Pair (NM_001128301)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57149
Forward Sequence:	ATCAGGGCAGATGGAAGACACC
Reverse Sequence:	TCAATCCTGGCTGTGCATTCATC
ACCN:	NM_001128301 , NM_001128301.1 , NM_001128301.2 , BC017039 , BC017039.1 , BF575136 , BG541469 , BI547471 , BI825882 , BM546065 , BX096767 , NM_001128301.3
UniProt ID:	O43325
Synonyms:	A211C6.1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)