

Product datasheet for **HP228309**

ATXN7L1 Human qPCR Primer Pair (NM_152749)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	222255
Forward Sequence:	GTGCCTGTAACCAGGTCGCAA
Reverse Sequence:	CCTGATGCTTCTCTATGGCTTGG
ACCN:	NM_152749 , NM_152749.1 , NM_152749.2 , NM_152749.3 , BC030616 , BC030616.2 , BC003517 , BC055009 , BP339600 , BQ227169 , BQ431139 , BQ436521
UniProt ID:	Q9ULK2
Synonyms:	ATXN7L4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)