

Product datasheet for **HP228319**

BPOZ (ABTB1) Human qPCR Primer Pair (NM_172027)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80325
Forward Sequence:	CCTGTCTCTGACGGCTTCAACA
Reverse Sequence:	TCTCGGAAGTGGTCATCCAGCA
ACCN:	<u>NM_172027</u> , <u>NM_172027.1</u> , <u>NM_172027.2</u> , <u>BC011858</u> , <u>BC011858.2</u> , <u>BF664138</u> , <u>BM671027</u> , <u>NM_172027.3</u>
UniProt ID:	<u>Q969K4</u>
Synonyms:	BPOZ; BTB3; BTBD21; EF1ABP; PP2259
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)