

Product datasheet for **HP228352**

CCR11 (ACKR4) Human qPCR Primer Pair (NM_178445)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51554
Forward Sequence:	GTCTCTGGAATGCAGTTTCTGGC
Reverse Sequence:	GGTATGCTCAGCAAGATGGCAG
ACCN:	<u>NM_178445</u> , <u>NM_178445.1</u> , <u>NM_178445.2</u> , <u>BC069438</u> , <u>BC069438.1</u> , <u>BC095501</u> , <u>BQ015927</u>
UniProt ID:	<u>Q9NPB9</u>
Synonyms:	CC-CKR-11; CCBP2; CCR-11; CCR10; CCR11; CCRL1; CCX-CKR; CCX CKR; CKR-11; PPR1; VSHK1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)