

Product datasheet for **HP228389**

BIN2 Human qPCR Primer Pair (NM_016293)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51411
Forward Sequence:	GGCTGCTATGTGACCATCTTCC
Reverse Sequence:	GCTTCTCCAGTTTGCTCATCACC
ACCN:	NM_016293 , NM_016293.1 , NM_016293.2 , NM_016293.3 , BC047686 , BC047686.1 , BG388638 , BX281612 , NM_016293.4
UniProt ID:	Q9UBW5
Synonyms:	BRAP-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)