

Product datasheet for **HP228420**

MOX1 (MEOX1) Human qPCR Primer Pair (NM_004527)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4222
Forward Sequence:	GAGATTGCGGTAAACCTGGACC
Reverse Sequence:	TCTGAACTGGAGAGGCTGTGG
ACCN:	<u>NM_004527</u> , <u>NM_004527.1</u> , <u>NM_004527.3</u> , <u>BC069474</u> , <u>BC069474.1</u> , <u>BC069506</u> , <u>BC069506.1</u> , <u>BE898972</u> , <u>BG677020</u> , <u>BQ003431</u> , <u>BQ028272</u> , <u>BX284090</u> , <u>NM_004527.4</u>
UniProt ID:	<u>P50221</u>
Synonyms:	KFS2; MOX1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)