

Product datasheet for **HP228475**

RPP14 Human qPCR Primer Pair (NM_007042)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11102
Forward Sequence:	CACAGACTGATGTGGCTACCTTC
Reverse Sequence:	GAGCTGAGATAAGTCCGTTGATC
ACCN:	NM_007042 , NM_007042.1 , NM_007042.2 , NM_007042.3 , NM_007042.4 , BC012017 , BC012017.1 , BC002441 , BC007342 , BP352444 , BP352480 , NM_007042.6
UniProt ID:	O95059
Synonyms:	P14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)