

## Product datasheet for **HP228484**

### Synaptopodin (SYNPO) Human qPCR Primer Pair (NM\_007286)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11346
Forward Sequence:	AGGAGGTGAGATGCAGCACACT
Reverse Sequence:	TAGGGTGTGGGCTGGATGTCA
ACCN:	<u><a href="#">NM_007286</a></u> , <u><a href="#">NM_007286.1</a></u> , <u><a href="#">NM_007286.2</a></u> , <u><a href="#">NM_007286.3</a></u> , <u><a href="#">NM_007286.4</a></u> , <u><a href="#">NM_007286.5</a></u> , <u><a href="#">BC042624</a></u> , <u><a href="#">BC113927</a></u> , <u><a href="#">BC114512</a></u> , <u><a href="#">BC142683</a></u> , <u><a href="#">BC146665</a></u> , <u><a href="#">BC152538</a></u> , <u><a href="#">BP237193</a></u> , <u><a href="#">NM_007286.6</a></u>
UniProt ID:	<u><a href="#">Q8N3V7</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)