

Product datasheet for **HP228548**

OPA3 Human qPCR Primer Pair (NM_025136)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80207
Forward Sequence:	CCGCTTGCCAACCGTATTAAGG
Reverse Sequence:	TCCTCGTTCAGCGCTTGATGA
ACCN:	<u>NM_025136</u> , <u>NM_025136.1</u> , <u>NM_025136.2</u> , <u>NM_025136.3</u> , <u>BC005059</u> , <u>BC005059.1</u> , <u>BC047316</u> , <u>BC064146</u> , <u>BU617706</u>
UniProt ID:	<u>Q9H6K4</u>
Synonyms:	MGA3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)