

Product datasheet for **HP228592**

CCDC102B Human qPCR Primer Pair (NM_024781)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79839
Forward Sequence:	TGAGACTGAAAGCAATAAATCTGCC
Reverse Sequence:	GACAAAGCCGACTCCAGTCTCT
ACCN:	<u>NM_024781, NM_024781.1, NM_024781.2, BC056269, BC126448, BC126450, BQ003970</u>
UniProt ID:	<u>Q68D86</u>
Synonyms:	ACY1L; AN; C18orf14; HsT1731
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)