

Product datasheet for **HP228608**

Angiomotin (AMOT) Human qPCR Primer Pair (NM_133265)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	154796
Forward Sequence:	TCCAATGCTGGATCAGGCTTGC
Reverse Sequence:	ATTCAGCACGGTAGTCTCCACC
ACCN:	<u>NM_133265</u> , <u>NM_133265.1</u> , <u>NM_133265.2</u> , <u>BC130294</u> , <u>BC039408</u> , <u>BC050581</u> , <u>BC094712</u> , <u>NM_133265.3</u>
UniProt ID:	<u>Q4VCS5</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)