

Product datasheet for **HP228677**

DLGAP1 Human qPCR Primer Pair (NM_004746)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9229
Forward Sequence:	ACCGTTAGGACCATCCAGAGCA
Reverse Sequence:	GTGCTACTCTGGGCTGTGATAG
ACCN:	<u>NM_004746</u> , <u>NM_004746.1</u> , <u>NM_004746.2</u> , <u>NM_004746.3</u> , <u>BC030096</u> , <u>BC040718</u> , <u>BC045785</u> , <u>BC070175</u> , <u>BC136453</u> , <u>BC136454</u> , <u>NM_004746.4</u>
UniProt ID:	<u>O14490</u>
Synonyms:	DAP-1; DAP-1-ALPHA; DAP-1-BETA; DAP1; DLGAP1A; DLGAP1B; GKAP; SAPAP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)