

## Product datasheet for **HP228694**

### KLC4 Human qPCR Primer Pair (NM\_201523)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	89953
Forward Sequence:	GCTGTGCTCTATGGCAAAGGG
Reverse Sequence:	GCCTCATACTTGCCTGGTTTTG
ACCN:	<u><a href="#">NM_201523</a></u> , <u><a href="#">NM_201523.1</a></u> , <u><a href="#">NM_201523.2</a></u> , <u><a href="#">BC080637</a></u> , <u><a href="#">BC080637.1</a></u> , <u><a href="#">BC012357</a></u> , <u><a href="#">BC103727</a></u> , <u><a href="#">BQ688778</a></u> , <u><a href="#">BX640717</a></u> , <u><a href="#">NM_201523.3</a></u>
UniProt ID:	<u><a href="#">Q9NSK0</a></u>
Synonyms:	bA387M24.3; KNSL8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)