

Product datasheet for **HP228743**

Epac1 (RAPGEF3) Human qPCR Primer Pair (NM_006105)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10411
Forward Sequence:	GTCATTTCTGCGTGTGGACAAG
Reverse Sequence:	CCACTTTGCCATGTTCTTCCAGC
ACCN:	<u>NM_006105</u> , <u>NM_006105.1</u> , <u>NM_006105.2</u> , <u>NM_006105.3</u> , <u>NM_006105.4</u> , <u>NM_006105.5</u> , <u>BC017728</u> , <u>BC040534</u> , <u>BC068477</u> , <u>BC092404</u> , <u>BQ007070</u>
UniProt ID:	<u>O95398</u>
Synonyms:	bcm910; CAMP-GEFI; EPAC; EPAC1; HSU79275
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)