

## Product datasheet for **HP228746**

### WDR62 Human qPCR Primer Pair (NM\_173636)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	284403
Forward Sequence:	CACTTGCTGGAGATTGACCACC
Reverse Sequence:	TCTGGCTCACACTCTTCCTCCA
ACCN:	<a href="#">NM_173636</a> , <a href="#">NM_173636.1</a> , <a href="#">NM_173636.2</a> , <a href="#">NM_173636.3</a> , <a href="#">NM_173636.4</a> , <a href="#">BC017261</a> , <a href="#">BC058939</a> , <a href="#">BC128417</a> , <a href="#">BC151839</a> , <a href="#">BQ880028</a> , <a href="#">BQ960310</a> , <a href="#">BU508347</a> , <a href="#">BX647726</a> , <a href="#">NM_173636.5</a>
UniProt ID:	<a href="#">O43379</a>
Synonyms:	C19orf14; MCPH2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)