

## Product datasheet for **HP228758**

### PCBP2 Human qPCR Primer Pair (NM\_005016)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5094
Forward Sequence:	ATCTGCTGCCAGCATTAGCCTG
Reverse Sequence:	GGTGGTGAACAGCAGAAAGGGA
ACCN:	<u><a href="#">NM_005016</a></u> , <u><a href="#">NM_005016.1</a></u> , <u><a href="#">NM_005016.2</a></u> , <u><a href="#">NM_005016.3</a></u> , <u><a href="#">NM_005016.4</a></u> , <u><a href="#">NM_005016.5</a></u> , <u><a href="#">BC001155</a></u> , <u><a href="#">BC001155.1</a></u> , <u><a href="#">BC071942</a></u> , <u><a href="#">BC071942.1</a></u> , <u><a href="#">BC107688</a></u> , <u><a href="#">BC035420</a></u> , <u><a href="#">BE748080</a></u> , <u><a href="#">BF509278</a></u> , <u><a href="#">BM764315</a></u> , <u><a href="#">BU844465</a></u> , <u><a href="#">NM_005016.6</a></u>
UniProt ID:	<u><a href="#">Q15366</a></u>
Synonyms:	hnRNP-E2; HNRNPE2; HNRPE2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)