

Product datasheet for **HP228777**

PTPN5 Human qPCR Primer Pair (NM_032781)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84867
Forward Sequence:	GCTTCATTGCCACCAGCATCTG
Reverse Sequence:	TGGTACTGCTCGCATGTCTGGA
ACCN:	NM_032781 , NM_032781.1 , NM_032781.2 , NM_032781.3 , BC064807 , BC064807.1 , BC039897 , BC046435 , NM_032781.4
UniProt ID:	P54829
Synonyms:	PTPSTEP; STEP; STEP61
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)