

Product datasheet for **HP228807**

ST3GAL3 Human qPCR Primer Pair (NM_174963)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6487
Forward Sequence:	GCCTGCTGAATTAGCCACCAAG
Reverse Sequence:	CCACTTGCGAAAGGAGTCATCC
ACCN:	<u>NM_174963</u> , <u>NM_174963.1</u> , <u>NM_174963.2</u> , <u>NM_174963.3</u> , <u>BC050380</u> , <u>BC050380.1</u> , <u>NM_174963.5</u>
UniProt ID:	<u>Q11203</u>
Synonyms:	DEE15; EIEE15; MRT12; SIAT6; ST3GALII; ST3Gal III; ST3GalIII; ST3N
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)