

Product datasheet for **HP228813**

NHP2 Human qPCR Primer Pair (NM_017838)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55651
Forward Sequence:	CCTCACGCGGAAGCTCTACAAA
Reverse Sequence:	CAGTGTGTCTCTGCCAAAACC
ACCN:	NM_017838 , NM_017838.1 , NM_017838.2 , NM_017838.3 , BC006387 , BC000009 , BE797477 , BG827491 , BM913344 , NM_017838.4
UniProt ID:	Q9NX24
Synonyms:	DKCB2; NHP2P; NOLA2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)