

## Product datasheet for **HP228923**

### ZMAT5 Human qPCR Primer Pair (NM\_019103)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55954
Forward Sequence:	GCCAAGAAGGTCTGGTACGACA
Reverse Sequence:	AGTTGGAGCCAAAGTCGCACTG
ACCN:	<u><a href="#">NM_019103</a></u> , <u><a href="#">NM_019103.1</a></u> , <u><a href="#">NM_019103.2</a></u> , <u><a href="#">BC009717</a></u> , <u><a href="#">BC009717.1</a></u> , <u><a href="#">BG285784</a></u> , <u><a href="#">BK005198</a></u> , <u><a href="#">BQ931527</a></u> , <u><a href="#">BU598307</a></u>
UniProt ID:	<u><a href="#">Q9UDW3</a></u>
Synonyms:	SNRNP20; U11/U12-20K; ZC3H19
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)