

Product datasheet for **HP228929**

ZNF638 Human qPCR Primer Pair (NM_014497)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27332
Forward Sequence:	CCTGAGATCCTCCCATCGAGAA
Reverse Sequence:	CGACGGAATCTGTGACTTGAGC
ACCN:	<u>NM_014497</u> , <u>NM_014497.1</u> , <u>NM_014497.2</u> , <u>NM_014497.3</u> , <u>NM_014497.4</u> , <u>BC010355</u> , <u>BC024000</u> , <u>BC064530</u> , <u>BC083513</u> , <u>BC143728</u> , <u>BM152722</u> , <u>BP348645</u> , <u>BX537425</u> , <u>NM_014497.5</u>
UniProt ID:	<u>Q14966</u>
Synonyms:	NP220; ZFML; Zfp638
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)