

Product datasheet for **HP228974**

THYN1 Human qPCR Primer Pair (NM_199298)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29087
Forward Sequence:	AGCCGCCTAGAGAAAGGTGTAG
Reverse Sequence:	CATGGCTCTAAGGAAGTTCCGAG
ACCN:	<u>NM_199298</u> , <u>NM_199298.1</u> , <u>BC006978</u> , <u>BC006978.1</u> , <u>BC093074</u> , <u>BQ640137</u> , <u>BQ959408</u> , <u>BU658940</u> , <u>BU685681</u> , <u>BU934600</u> , <u>NM_199298.2</u>
UniProt ID:	<u>Q9P016</u>
Synonyms:	HSPC144; MDS012; MY105; THY28; THY28KD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)