

Product datasheet for **HP228979**

PARP6 Human qPCR Primer Pair (NM_020214)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56965
Forward Sequence:	GCCTATGGCAAAGGCATCTACC
Reverse Sequence:	TCTCTGGACCAGCTCATCCTTG
ACCN:	<u>NM_020214</u> , <u>NM_020214.1</u> , <u>NM_020214.2</u> , <u>NM_020214.3</u> , <u>BC110902</u> , <u>BC110902.1</u> , <u>BC015176</u> , <u>BC015314</u> , <u>BC026955</u> , <u>BC119688</u> , <u>BC119689</u> , <u>NM_020214.4</u>
UniProt ID:	<u>Q2NL67</u>
Synonyms:	ARTD17; PARP-6-B1; PARP-6-C; pART17
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)