

## Product datasheet for **HP228985**

### AP3M1 Human qPCR Primer Pair (NM\_207012)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26985
Forward Sequence:	CATACAGACCGAAGTGCCACCT
Reverse Sequence:	CCTTAATTGCAGCCTCTGAACAC
ACCN:	<a href="#">NM_207012</a> , <a href="#">NM_207012.1</a> , <a href="#">NM_207012.2</a> , <a href="#">NM_207012.3</a> , <a href="#">BC067127</a> , <a href="#">BC026232</a> , <a href="#">BE706537</a> , <a href="#">BM463185</a> , <a href="#">BP230128</a> , <a href="#">BU172918</a> , <a href="#">BU430175</a> , <a href="#">BU679701</a> , <a href="#">NM_207012.4</a>
UniProt ID:	<a href="#">Q9Y2T2</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)