

Product datasheet for **HP228995**

DDX4 Human qPCR Primer Pair (NM_001136034)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54514
Forward Sequence:	GTGTCTGGACATGATGCCACCAC
Reverse Sequence:	GCAAGCCATCAAATCTCGTCCTG
ACCN:	<u>NM_001136034</u> , <u>NM_001136034.1</u> , <u>BC088362</u> , <u>BC088362.1</u> , <u>BC047455</u> , <u>BC030638</u> , <u>BI559479</u>
UniProt ID:	<u>Q9NQI0</u>
Synonyms:	DEAD (Asp-Glu-Ala-Asp) box polypeptide 4; DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide 4; MGC111074; OTTHUMP00000122546; VASA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)